

○

**READY. SET.
KEEP YOUR
LAB ON THE
GO.**

Streamline media preparation with
Thermo Scientific Bagged Enrichment Media

Thermo Scientific Bagged Enrichment Media

Prepare to be more productive with our range of convenient, ready-to-hydrate solutions

Increase productivity by decreasing preparation time with Thermo Scientific™ Bagged Enrichment Media. Containing pre-weighed, gamma-irradiated enrichment media, our scalable solutions enable you to respond rapidly to meet fluctuating testing demands and focus resources on higher value activities—keeping your laboratory on the go.



Fast

Prepare high quality enrichment media in a matter of minutes.



Flexible

Work with confidence, knowing that your enrichment media will be ready to go when you are.



Hassle-free

Eliminate weighing, autoclaving and bottlenecks from your workflow.



More formats to fit your needs

Choose from FitBag, QuickBag or Dry-Bags Media

Designed to meet your sample preparation and testing volume needs, Thermo Scientific Bagged Enrichment Media comes in three easy-to-use formats, keeping medium preparation convenient and cost-effective.



FitBag Media

Save space and time during medium preparation with easy-to-use Thermo Scientific™ FitBag™ Media.

Just fill with water*, break the dividing seal, and use as needed.

QuickBag Media

Reduce the risk of contamination during enrichment medium preparation with Thermo Scientific™ QuickBag™ Media.

To use, simply break the dividing seal, mix, and start dispensing.

Dry-Bags Media

Make large-scale enrichment medium preparation more efficient with Thermo Scientific™ Dry-Bags™ Media.

Add water* to the pre-weighed powder, mix, and dispense.

Preparation time

< 10 minutes

< 5 minutes

< 30 minutes

Available volumes

2.7 L, 4.5 L, 9 L

2.7 L, 4.5 L

20 L

Features include

Flexibility to fill when convenient and mix as needed**

Pre-filled water compartment***


Large volume format for high test throughput***

FitBag Media can be efficiently stored, requiring less storage space than ready-to-use media in bottles and bags.

5
Media bottles



=



50
Un-filled FitBag
Media Bags

* Use a Thermo Scientific Peristaltic pump and Thermo Scientific 0.2µ Filter to fill FitBag and Dry-Bags Media with deionized water.

** FitBag Media bags can be filled up to one week before mixing. Once hydrated FitBag Media should be used within 72 hours.

*** Once hydrated FitBag and Dry-Bags Media should be used within 72 hours.

Complete solutions. Complete convenience.

Ensure efficiency across your entire workflow with Thermo Scientific sample preparation and enrichment solutions



Peristaltic Pumps

Fill your FitBag and Dry-Bags Media the quick and easy way with Thermo Scientific™ Peristaltic Pumps.



DiluFlux Pro Automated Gravimetric Dilutors

Dilute samples quickly and accurately at the touch of a button with the Thermo Scientific™ DiluFlux Pro™ Dilutors.



Homogenizer Laboratory Blender

Achieve powerful, quiet sample homogenization in under a minute using the Thermo Scientific™ Homogenizer Laboratory Blender.

Prepare to be more productive with Thermo Scientific sample preparation and enrichment solutions. Our comprehensive range of high-quality products are designed to be fast, flexible and freeing – enabling you to meet testing demands with convenience and confidence.

To learn more, please contact your local sales representative or visit thermofisher.com/enrich

© 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. FitBag and QuickBag are trademarks of Interbiotech S.r.l.

Contact Information:

microbiology@thermofisher.com
USA +1 800 255 6730
International +44 (0) 1256 841144

993-002
LT2340A
May 2017

ThermoFisher
SCIENTIFIC