

Case study

Easily label increased sample volumes



A biological sample testing laboratory can label tubes much faster and more reliably with the BradyPrinter i3300 and B-492 FreezerBondz labels.

Challenge

How to identify thousands of blood samples efficiently

A biological sample testing laboratory needed to check blood samples for COVID-19 anti-bodies during the 2020 - 2021 pandemic.

To make this happen, huge quantities of blood samples had to be serologically tested and stored in tubes in freezers and then for up to 10 years in liquid nitrogen as a back-up. With huge and growing numbers of infected people, the process had to be fast and efficient.

One of the challenges the laboratory needed to address quickly was sample labelling. Their current label printer was failing them and absorbed too much time because it had to be calibrated manually.

Solution

Auto-calibration and software integration

Together with a specialised distributor, Brady proposed the BradyPrinter i3300 Industrial Label Printer together with B-492 labels. The printer recognises the label material on insertion, immediately auto-calibrates and is ready to deliver perfect prints within 20 seconds. This significantly reduced time spent by scientists on label printing, time much better spent on testing blood samples in the fight against the pandemic.

Additionally, respecting anti-viral social distancing measures, a call was set-up to integrate the printer with the customer's existing software systems. This resulted in automated label printing within an existing process scientists were already familiar with.



On top of this, the B-492 'FreezerBondz' label that could now be printed almost automatically, with the BradyPrinter i3300 isn't just any label. It perfectly fits tubes, can resist common laboratory chemicals, and the label stays legible and attached under liquid nitrogen storage temperatures at -196°C .

Result

Fast, right first-time label printing

Scientists in the biological sample testing laboratory can now quickly print and apply a reliable label on a blood sample tube. With printer auto-calibration, they no longer need to spend time on printer settings and re-printing labels.

Visit eu.fishersci.com for more information

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2024 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

23094

 **fisher scientific**
part of Thermo Fisher Scientific